

Proven™ Revision Knee System

Product Specifications



Proven™ Revision Knee System

Product Specifications

StelKast's patented Proven Revision Knee System features unsurpassed stem offset positioning.

Revision Femur

- Adjustable A-P Stem Location aids in precise component fit and sizing.
- Multiple Valgus Angles – 5 & 7 Degrees Standard
- Mechanically Locked Distal & Posterior augmentation

Revision Tibial Tray

- Adjustable Tibial Stem Offsets – 0, 2 & 4mm ensures proper fit and bone coverage
- Proximal Offset places Stem Offset Close to Tibial Plateau
- Mechanically Locked Block Augmentation
- Tri-Slot Titanium Cementless Stems
- Polished Cobalt-Chrome Cemented Stems

Proven Revision Knee System		
Catalog#	Description	Offerings
SC2957	Revision Femur	Left – Sizes 1 thru 5
SC2958		Right – Sizes 1 thru 5
SC2193	Valgus Bushing	5 & 7 Degrees Standard (0 & 3 Degrees optional)
SC2316	Distal Femoral Augment	Sizes 1 thru 5, 5 & 10mm Thickness
SC2317	Posterior Femoral Augment	Sizes 1 thru 5, 5 & 10mm Thickness
SC3029	Femoral & Tibial Stem	Cemented – 10-16mm Diameter, 40, 80 & 120mm Length
SC2196		Cementless – 10-24mm Diameter, 80 & 120mm Length
SC2241	Revision Tibial Tray	1T1F, 1T2F, 2T2F, 2T3F, 3T2F, 3T3F, 3T4F, 4T3F, 4T4F, 4T5F, 5T4F, 5T5F
SC2240	Revision Tibial Bushing	0, 2 & 4mm Offset
SC2287	½ Tibial Block (w/screws)	Sizes 1 thru 5, 5 & 10mm Thickness
SC2265	Revision Insert	Sizes 1 thru 5, 10, 12, 14, 17, 20, 23, 26, & 30mm Thickness
SC2291	Modular Tibial Tray	1T1F, 1T2F, 2T2F, 2T3F, 3T2F, 3T3F, 3T4F, 4T3F, 4T4F, 4T5F, 5T4F, 5T5F



200 Hidden Valley Road | McMurray, PA 15317