

GENFlex₂

Brochure / Product Specifications



GENFlex₂ Brochure / Product Specifications

- Implant designs allow the freedom and range of motion expected by today's active patients
- 10 femoral sizes address diverse patient population
- Implant fit optimized using StelKast's extensive bone database combined with advanced simulation and virtual implantation
- Insert thicknesses in 1mm increments for precision gap balancing
- EXp Blended Vitamin E polyethylene options for tibial inserts and patella components
- Anterior flange angle enables downsizing of femur for gap balancing while reducing the incidence of femoral notching
- Comprehensive total knee system allows efficient flow from primary to revision knee conversion
- Multiple instrumentation options accommodate various surgical techniques

Catalog #	Description	Offerings
SC3749	Posterior Stabilized (PS) Femur, Left	Sizes 1*, 1.5, 2, 2.5, 3, 3.5, 4, 4.5, 5, 6
SC3750	Posterior Stabilized (PS) Femur, Right	Sizes 1*, 1.5, 2, 2.5, 3, 3.5, 4, 4.5, 5, 6
SC3751	Cruciate Retaining (CR) Femur, Left	Sizes 1*, 1.5, 2, 2.5, 3, 3.5, 4, 4.5, 5, 6
SC3752	Cruciate Retaining (CR) Femur, Right	Sizes 1*, 1.5, 2, 2.5, 3, 3.5, 4, 4.5, 5, 6
SC3753	PS Insert	Sizes 1* thru 6 10, 11, 12, 13, 15, 17, 20, 23*, 26*, 30*mm
SC3754	CR Insert	Sizes 1* thru 6 10, 11, 12, 13, 15, 17, 20mm
SC3802	EXp PS Insert	Sizes 1* thru 6 10, 11, 12, 13, 15, 17, 20, 23*, 26*, 30*mm
SC3803	EXp CR Insert	Sizes 1* thru 6 10, 11, 12, 13, 15, 17, 20mm
SC3915	CR Ultracongruent Insert	Sizes 1* thru 6 10, 11, 12, 13, 15, 17, 20, 23*, 26*, 30*mm
SC3916	EXp CR Ultracongruent Insert	Sizes 1* thru 6 10, 11, 12, 13, 15, 17, 20, 23*, 26*, 30*mm
SC1584	Tibial Tray	1T1F, 1T2F, 2T1F, 2T2F, 2T3F, 3T2F, 3T3F, 3T4F, 4T3F, 4T4F, 4T5F, 5T4F, 5T5F, 5T6F, 6T5F, 6T6F
SC1591	3 Peg Patella	26, 29, 32, 35, 38, 41mm
SC1592	1 Peg Patella	26, 29, 32, 35, 38, 41mm
SC3454	EXp 3 Peg Patella	26, 29, 32, 35, 38, 41mm
SC3426	EXp 1 Peg Patella	26, 29, 32, 35, 38, 41mm

* Indicates Special Order.

www.stelkast.com



200 Hidden Valley Road | McMurray, PA 15317
phone: (724) 941-6368 | fax: (724) 941-5987